

**URGENT FIELD SAFETY NOTICE****Potential for False Positive Results Using ORTHO BioVue® System (Rh/K Cassette), Lot RHP020F****Date Issued****xx July 2017****Issue**

As part of a Field Safety Corrective Action, this Urgent Field Safety Notice is in reference to the following products:

Product Name	Product Code	Lot Number	Expiry
Blood Grouping Reagents Anti-C (Anti-RH2) (Monoclonal) Anti-E (Anti-RH3) (Monoclonal) Anti- $\bar{c}$ (Anti-RH4) (Monoclonal) Anti-e (Anti-RH5) (Monoclonal) Anti-K ORTHO BioVue® System (Rh/K Cassette)*	707280 (400 Cassettes) 707250 (100 Cassettes)	RHP020F	18-OCT-17

\*ORTHO BioVue System (Rh/K Cassette) is a qualitative test for recognition of the C (RH2), E (RH3),  $\bar{c}$  (RH4), e (RH5) and K (K1) antigens on human red blood cells

Ortho has confirmed the potential for unexpected false positive 3+ or 4+ results for either anti- $\bar{c}$  (column 3) or anti-e (column 4) in a limited number of ORTHO BioVue® System (Rh/K Cassettes), Lot RHP020F only. The unexpected results are caused by a dried column that is located in either the anti- $\bar{c}$  (column 3) or anti-e (column 4). If the cassette is affected, the dried column is located in either the anti- $\bar{c}$  (column 3) or anti-e (column 4), but not both.

Ortho determined that up to 0.6% of 405,000 total cassettes produced for this lot only may be affected. If an affected cassette is used for testing:

- On ORTHO AutoVue® *Innova* or *Ultra* Systems, unexpected false positive results may occur without operator notification.
- On ORTHO VISION® and ORTHO VISION® Max Analyzers, the system automatically detects dried columns and rejects affected cassettes prior to using them for testing.
- In manual BioVue cassette testing, unexpected false positive results may occur.

Please refer to Questions and Answers on Page 3 for additional information.

**Required Actions**

- Immediately discontinue using and discard your remaining inventory of ORTHO BioVue System (Rh/K Cassette), Lot RHP020F. **NOTE:** If you have no alternative lot in your inventory, you may use the affected lot until you receive your replacement order if:
  1. Using Ortho VISION or ORTHO VISION Max Analyzers or
  2. You manually inspect cassettes from the affected lot for dried columns prior to use if performing:
    - a. manual BioVue cassette testing, or
    - b. testing on ORTHO AutoVue® *Innova* or *Ultra* Systems.

---

**Required  
Actions,  
continued**

- Complete and return the Confirmation of Receipt form no later than **xx July 2017**. Upon receiving your Confirmation of Receipt form, Ortho will credit your account or send a replacement order for discarded product.
  - Forward this notification if you have provided this product outside of your facility.
- 

**Impact to  
Results**

Use of cassettes with either a dried anti- $\bar{c}$  column or dried anti-e column for testing could potentially mistype patients who are antigen negative as being  $\bar{c}$  or e antigen positive and erroneously suitable to receive  $\bar{c}$  or e positive blood. If the cassette is affected, the dried column is located in either the anti- $\bar{c}$  (column 3) or anti-e (column 4), not both.

Please review previous results generated from testing with this product lot and consult with your Medical Director to determine whether any follow-up activities are required for samples that produced 3+/4+ anti- $\bar{c}$  or anti-e results only for patients who have no previous clinical history.

If visual inspection was performed prior to use, and dried columns were identified and rejected for use, previous results utilizing this lot should not be affected. As a reminder, the Precautions Section of the Instructions for Use currently states:

*Do not use cassettes that appear damaged (i.e., break in foil seal or break, crack or bubble in the column) or exhibit drying (i.e., liquid level is at or below the top of the glass beads) or exhibit discoloration (due to bacterial contamination which can cause false reactions).*

---

**Investigation**

Ortho became aware of this quality issue from customers who contacted us after either finding dried columns when manually checking cassettes prior to use or after obtaining unexpected false positive results when using their ORTHO AutoVue *Innova* or *Ultra* Systems.

Ortho's root cause investigation in coordination with our plastic cassette supplier has traced this issue to a defect in the plastic molding of a limited number of cassettes produced at the supplier's facilities. Only lot RHP020F is affected by this issue.

---

**Contact  
Information**

We apologize for the inconvenience this will cause your laboratory. If you have further questions, please contact Ortho Care™ Technical Solutions Center at **insert appropriate number / insert signatory if required**

---

Enclosure: Confirmation of Receipt Form

## Questions and Answers

### 1. Is there any impact to previously reported results if I used cassettes from the affected lot?

A dried column could cause a false positive or erroneous result. If the cassettes from Lot RHP020F were inspected prior to use and the result was as expected, no further action is required. If you suspect that an erroneous result occurred (i.e., 3+/4+ anti-c̄ or anti-e reaction only for patients who have no previous clinical history), please consult with your Medical Director and report the occurrence to our Ortho Care™ Technical Solutions Center at *insert appropriate number*.

### 2. Are any other products affected by this issue?

Ortho has determined that this issue is limited to 0.6% of the cassettes manufactured from ORTHO BioVue System (Rh/K Cassette), Lot RHP020F only. The affected cassettes contain a plastic molding defect that occurred during a single, identified manufacturing event that occurred at the supplier of our plastic cassettes.

### 3. Will this issue be detected if I am using an ORTHO VISION Analyzer or ORTHO VISION Max Analyzer or BioVue manual testing?

As part of their design, ORTHO VISION and VISION Max Analyzers perform a quality pre-check and reject cassettes identified with dried columns (Error Code CIMS21) prior to use. Unlike ORTHO VISION and ORTHO VISION Max Analyzers, ORTHO AutoVue *Innova* and *Ultra* Systems do not have the functionality to detect dried columns prior to cassettes being used for testing.

The ORTHO BioVue System (Rh/K Cassette) IFU cautions against using cassettes for testing that exhibit drying (i.e., liquid level is at or below the top of the glass beads).

### 4. When can I expect my replacement order?

For your convenience, Ortho will either process your product replacement or credit your account upon receipt of your Confirmation of Receipt Form. Please choose your replacement or credit preference on the Confirmation of Receipt Form included with this communication.

### 5. I do not have an alternative lot in my inventory. May I continue to use the affected lot until I receive my replacement order?

Yes. If you have no alternative lot in your inventory, you may use the affected lot until you receive your replacement order if:

1. Using Ortho VISION or ORTHO VISION Max Analyzers, or
2. You manually inspect cassettes from the affected lot for dried columns prior to use if performing:
  - a. manual BioVue cassette testing, or
  - b. testing on ORTHO AutoVue® *Innova* or *Ultra* Systems.