


Attachment: Instruction for manual setting of sample probe washing on cobas c 503

Step 1: Select "Menu"



Step 2: Select "Special Wash"



Step 3: Select "Sample Probe" for c 503 module



Step 4: In section "User rules" select empty line and select button on the right and "Edit"

User Rules			
Select	Victim	Washing Method	Material
<input type="checkbox"/> 1			
<input type="checkbox"/> 2			

Step 5: Search and select IGG2C test, ACN # 20741

Select

All Select

Search Criteria Test ACN

Test IGG2C

ACN Test

20741 IGG2C

Step 6: Select Washing Method "Wash Solution"

Assign Wash Rule to Sample Probe

Rule No.: 1

Module: c 503

Victim: IGG2C

Washing Method: Wash Solution

Material: Samples QC Calibrator

Step 7: Assign wash rule: Samples and QC

Assign Wash Rule to Sample Probe

Rule No.: 1

Module: c 503

Victim: IGG2C

Washing Method: Wash Solution

Material: Samples QC Calibrator

Step 8: Save