

Siemens Healthcare Diagnostics GmbH, SHS EMEA CEET QT, Siemensstrasse 90,  
1210 Vienna

Name M.A. Roland Ertl  
 Department SHS EMEA CEET QT  
 Mobile +43 (664) 8011738274  
 E-mail roland.re.ertl@siemens-healthineers.com  
 Date January 09, 2024  
 Document Ref# ACHC24-01.A.OUS

**Urgent Field Safety Notice:**

**Atellica® CH Analyzer**  
**Atellica® CI Analyzer**

**Potential for Negative Bias with Atellica CH Immunoglobulin M\_2 (IgM\_2) Reagent**

– To whom it may concern,

Our records indicate that your facility may have received the following product:

**Table 1. Atellica CH and Atellica CI Affected Product**

Assay	Test Code	Siemens Material Number (SMN)	Unique Device Identification (UDI)	Lot Number
Immunoglobulin M_2	IgM_2	11097620	00630414595627	221764 and above

**Reason for Correction**

The purpose of this communication is to inform you of an issue with the product indicated in Table 1 above and provide instructions on actions that your laboratory must take. Siemens Healthineers has confirmed the potential for a negative bias with quality control (QC) and patient sample results when using the Atellica CH Immunoglobulin M\_2 (IgM\_2) reagent. The negative bias was observed after the IgM\_2 reagent was stored onboard the analyzer regardless of whether the reagent wells were punctured or unpunctured. Unopened reagents stored refrigerated at 2 - 8 °C are unaffected. See Tables 2 and 3 in the Additional Information section for observed results.

This correction is applicable to all future lots until further notice. Siemens Healthineers is currently investigating the root cause of this issue.

**Risk to Health**

When this issue occurs, there is a potential for erroneously depressed IgM patient results. This is not expected to lead to a significant effect on assessment of IgM results in the context of the assay intended use.

## Actions to be Taken by the Customer

Siemens Healthineers recommends batch testing samples for Atellica CH IgM\_2 as follows:

1. Remove and discard any Atellica CH IgM\_2 reagent packs onboard the analyzer.
2. Load a **single** fresh Atellica CH IgM\_2 reagent pack onto the analyzer.
3. Perform a **Lot calibration** and process Quality Control (QC).
4. Immediately process a batch of patient samples and conclude with a repeat run of QC.

Patient results should not be reported until the QC performed at the end of the batch run has been assessed.

- If the QC results are within the established range, patient results can be reported.
  - If the QC results are not within the established range, do not report patient results and repeat steps 1 – 4 above.
5. Remove and discard the Atellica CH IgM\_2 reagent pack at the end of the batch run.
- Siemens Healthineers does not recommend using the ADVIA IgM\_2 reagent on the Atellica CH or Atellica CI Analyzers.
  - Please review this letter with your Medical Director to determine the appropriate course of action, including for any previously generated results, if applicable.
  - Complete and return the Field Correction Effectiveness Check Form attached to this letter within 30 days.
  - If you have received any complaints of illness or adverse events associated with the product listed in Table 1, immediately contact your local Siemens Healthineers Customer Care Center or your local Siemens Healthineers technical support representative.

Please retain this letter with your laboratory records and forward this letter to those who may have received this product.

We apologize for the inconvenience this situation may cause. If you have any questions, please contact your Siemens Healthineers Customer Care Center or your local Siemens Healthineers technical support representative.

Sincerely yours,

Siemens Healthcare Diagnostics GmbH

Signature:



*Electronically signed by: Roland Ertl  
Reason: I am approving this document  
Date: Jan 8, 2024 10:37 GMT+1*

Email: roland.re.ertl@siemens-healthineers.com

i.A. Roland Ertl, MA  
Quality Management CEECA

Signature:



*Electronically signed by: Carina Marie  
Viehboeck  
Reason: I have reviewed this document  
Date: Jan 8, 2024 10:44 GMT+1*

Email: carina-marie.viehboeck@siemens-healthineers.com

i.A. DI<sup>in</sup> Carina Viehböck  
Product Manager CEECA

## Annex I: Additional Information

Siemens Healthcare Diagnostics GmbH  
Management: Joachim Bogner, Stefan Scheidler,  
Sonja Wehsely

Siemensstrasse 90  
1210 Vienna  
Austria

siemens-healthineers.com/at

## Annex I: Additional Information

**Table 2. Atellica CH IgM\_2 Quality Control Results Across Four Hours\***

Time (hour)	QC Level 1		QC Level 2		QC Level 3	
	mg/dL (g/L)		mg/dL (g/L)		mg/dL (g/L)	
	Rep 1	Rep 2	Rep 1	Rep 2	Rep 1	Rep 2
0 (after calibration)	47.0 (0.47)	45.5 (0.46)	76.3 (0.76)	72.2 (0.72)	90.4 (0.90)	88.9 (0.89)
1	41.6 (0.42)	38.9 (0.39)	72.2 (0.72)	70.1 (0.70)	79.8 (0.80)	84.7 (0.85)
2	38.4 (0.38)	37.1 (0.37)	65.4 (0.65)	64.7 (0.65)	77.5 (0.78)	78.5 (0.79)
3	34.5 (0.35)	33.9 (0.34)	63.3 (0.63)	64.7 (0.65)	75.5 (0.76)	75.4 (0.75)
4	32.2 (0.32)	33.5 (0.34)	63.1 (0.63)	63.1 (0.63)	73.2 (0.73)	75.9 (0.76)

\*Quality Control (QC) materials are representative of patient samples. Bio-Rad Multiquel Assayed Control Lot 45960 was used for testing.

**Table 3. Atellica CH IgM\_2 Patient %Bias Results at 24 Hours\*\***

Range of Patient Serum Sample IgM_2 Values	%Bias Range	Average %Bias
23 – 50 mg/dL (0.23 - 0.50 g/L)	[-14.9 to -26.5%]	-20.7%
51 – 100 mg/dL (0.51 - 1.00 g/L)	[-6.7 to -15.4%]	-11.0%
101 – 200 mg/dL (1.01 - 2.00 g/L)	[-4.1% to -6.3%]	-5.2%
201 – 330 mg/dL (2.01 - 3.30 g/L)	[-2.9% to -4.7%]	-3.8%

\*\*Initial result was obtained immediately after calibration of a freshly loaded reagent pack. A repeat result was obtained 24 hours after calibration from the same reagent well.