



URGENT FIELD SAFETY NOTICE

RECALL of CRP FS

Deviation in recovery in lower concentration range

Date 2025-10-02

Product / Lot CRP FS Cat. No.: 1 7002 ...
Lot: 35872/ 50308370/1, 50309017/8, 50312566,
Kit Lot: 60166351, 60166375, 60166095, 60166128,
60166320, 60165831, 60166238, 60166105,
60166148, 60166196, 60167664, 60166543,
60165674, 60165433, 60165899

Lot: 36176/ 50310306/7, 50311461/2, 50311720/1,
50310374
Kit Lot: 60168245, 60166682, 60167113, 60166458,
60167194, 60166381, 60166761, 60167251,
60167371, 60166725, 60167234, 60167518,
60167658, 60168139, 60167574, 60167496,
60166829, 60167320 60167124, 60167168,
60167123

Explanation Based on the analysis of a customer complaint the a.m. lots of CRP FS have been reviewed. The internal inspection confirms incorrect recovery of low concentrated samples (<10 mg/L), which leads to false low results. Current findings suggest that a specific raw material could be the cause.

Impact on patient results CRP FS measurements of patient samples lead to deviating results.

Measures The a.m. reagent lots are withdrawn from the market.

Please immediately inform all involved users about this recall. Immediately stop further measurements and replace the reagent lots at the customer's site. Existing stock kits must be disposed of. DiaSys offers replacement for the affected kits.

Determinations performed using a.m. lots of CRP FS shall be assessed with care and repeated if implausible results occur.

DiaSys has announced this field safety notice to the relevant authorities of the European Union. Customers outside the EU are asked to handle necessary announcements to authorities in their countries.



URGENT FIELD SAFETY NOTICE

Under current regulations, we are obliged to provide a complete chain of evidence of all corrective measures for our products. For this reason, we ask you to fill in and sign the attached confirmation of receipt as proof that you have received and communicated this information to all concerned customers. Please return it by fax or as scan until 2025-10-17.

Please accept our sincere apologies for the inconvenience caused. In case you have any questions, please do not hesitate to contact us.

Best regards,



Tim Fleissner

Strategic Product Management Reagents